NOVUS® Valves and Kits

Sterile For Single Use Only

STERILE

EO

0123

ENGLISH

NL850-9010	NL850-9020	NL850-9030	NL850-9011	NL850-9021	NL850-9031
NL850-9610	NL850-9620	NL850-9630	NL850-9311	NL850-9321	NL850-9331


Figure 1. **NOVUS® Mini Valve**


Figure 2. **NOVUS® Valve**

Description

The NOVUS® and NOVUS® Mini Valves are silicone elastomer valves used for the treatment of hydrocephalic patients when shunting cerebrospinal fluid (CSF) from the ventricles of the brain. The valves incorporate an internal low, medium or high pressure valve assembly which is mounted distal to an integral pumping reservoir for proximal control of CFSF flow. The devices also provide one-way flow control. An integral occluder is incorporated on the proximal end of the valves, while an integral Physiological Flow Device is mounted on the distal end and also serves as an occluder. The occluder and Physiological Flow Device are designed to manually control the flushing or pumping of CSF proximally or distally. The devices are fitted with a polypropylene needle guard, designed to prevent inadvertent needle puncture during percutaneous injection. The valves contain an integral plastic connector on each end to simplify the assembly procedure. The valves contain no metal parts to interfere with computerized axial tomography. Due to the internal reservoir capacity, no separate reservoir is usually required (Figures 1 and 2).

The Integra NeuroSciences Hydrocephalic Shunt Kits are designed to deliver CSF from the ventricles of the

U.S. Patent Nos. 4,781,674, 4,904,236, 3,769,982
4,364,395 and 4,464,168

NOVUS and NOVUS Mini valves are registered trademarks of Integra NeuroSciences.

102532-701-03

brain to the peritoneal cavity. Each kit contains either a NOVUS or NOVUS Mini Valve, a Pudenz Ventricular Catheter and an integral, open-ended kink resistant catheter attached to the outlet of the valve.

The kits are available in low, medium and high pressure ranges. All valves and catheters contain barium sulfate for x-ray detectability. No metal parts are used anywhere in the system, to eliminate any possibility of interference with CT scanning or MR imaging.


NOVUS Models

The NOVUS Models have a contoured base which facilitates choice of placement remote from the burr hole. A silicone elastomer T-valve within the assembly is designed to provide shunt pressure control and to help prevent backward flow of CSF into the ventricles of the brain.

Three valve closing pressure ranges are available to provide the surgeon with a choice in meeting individual patient needs. Each valve is individually tested at the time of manufacture for conformance with labeled closing pressure performance characteristics. The closing pressure ranges are identified by a tantalum-impregnated silicone elastomer dot code on the valve dome.

Valve performance characteristics are reported in accordance with a shunt valve test method contained in the American Society for Testing and Materials F647-86, "Standard Practice for Evaluating and Specifying Implantable Shunt Assemblies for Neurosurgical Application" and by valve closing pressure for the convenience of performance verification prior to implantation by the surgeon. For a description of the preimplant check for valve patency and valve performance, please refer to Instructions for Use. Due to characteristics of silicone materials, some variation in pressure performance may occur.

In order to allow the valve assembly to control shunt pressure, Integra NeuroSciences recommends the use of a low-pressure distal catheter (closing pressure 54mm H₂O or less) such as the Integra NeuroSciences Peritoneal Reflux Control Catheter (Catalog No. NL850-1375) or Peritoneal Open-Ended

Ventricular Catheter Guide

Catheter (Catalog No. NL850-1376).
See Instructions for Use.

The valve flushing chamber may be penetrated by a 25-gauge or smaller needle. Such penetration permits removal of CSF samples, as well as the injection of fluid materials in either the proximal or distal direction. The flushing chamber is marked with an arrow to indicate the CSF flow direction, which facilitates proper valve placement.

Physiological Flow Device

All models possess a Physiological Flow Device which also serves as distal occluder. This device is designed to help prevent the excessive drainage of cerebrospinal fluid which would otherwise be caused by the siphoning effect created by the elevation of the ventricular catheter with respect to the distal catheter (i.e. when the patient sits, stands or is held erect). This siphoning effect is minimized by the physiological flow device, which closes under normal physiological negative pressure, yet will reopen to allow the flow of CSF to resume before intraventricular pressure becomes excessive. The Physiological Flow Device is designed to restrict flow when the shunt system is under substantial negative hydrostatic pressure. It is not designed to compensate for lowered intracranial pressure brought on by other than large negative hydrostatic pressure changes.

Indications

The **NOVUS** and **NOVUS Mini** Valves utilized in the treatment of hydrocephalic patients are components in systems designed to shunt cerebrospinal fluid from the lateral ventricles of the brain into either the right atrium of the heart or the peritoneum.

The **NOVUS Mini** valve can be used in (but is not restricted to) situations

where skin erosion may be a problem, as with older patients.

Valves with a Physiological Flow Device are intended to reduce the hazard of negative intraventricular pressure (with respect to atmospheric pressure) when the patient is sitting, standing or semirecumbent.

Contraindications

Ventriculoatrial or ventriculoperitoneal shunting systems should not be used in the presence of known or suspected infections along the course of the shunt (meningitis, ventriculitis, skin infections, bacteremia, septicemia, or peritonitis). It is advisable to avoid shunting procedures if infection is present anywhere in the body.

The ventriculoatrial method of shunting is contraindicated for patients with congenital heart disease or other anomalies of the cardiopulmonary system.

Instructions for Use

Valve Patency and Closing Test

Note: There is no basic, concise test which can be performed under operating room conditions to measure specific valve pressure. The following test is provided only to confirm patency prior to implantation.

1. Connect a length of sterile tubing (approximately 25cm) to the proximal connector. No tubing should be connected to the outlet tube.
2. Orient the valve horizontally on the test table.
3. Using a syringe, gently fill the tube. Use sterile deaerated water to help eliminate air bubbles in tubing and valve.
4. After flow through tubing and valve has been established, examine tubing and valve to insure bubbles have been eliminated. **Do not flush the valve with a syringe.** The amount of pressure created by syringe flushing will temporarily deform the silicone elastomer valve, which may cause abnormally low pressure test results.
5. Elevate the tubing vertically above the valve. A noticeable slowing of flow

should occur when the water level approaches or falls within the indicated pressure range. The pressure is that recorded by measuring the distance from the base of the valve to the top of the water column in the proximal tube extension. Due to characteristics of silicone materials, some variation in pressure performance may occur.

A radiopaque ventricular catheter guide is provided to help direct the catheter out of the trephine opening to minimize kinking.

The **Pudenz Ventricular Catheter** is manufactured from silicone elastomer with a barium sulfate impregnated stripe to provide radiopacity throughout the entire length of the catheter. The catheter has multiple small perforations near the tip to facilitate drainage of CSF and to help reduce clogging of the catheter lumen. The extreme tip of the catheter is closed to help prevent particulate matter from entering the tubing during placement. The tip of the catheter contains a radiopaque, tantalum-impregnated silicone elastomer plug to aid in determining location during and subsequent to implantation. The catheter is provided with a stylet to support the catheter during implantation.

The peritoneal catheter which is integral to the outlet of the valve for some models is made of a high durometer (firm) silicone elastomer which helps resist kinking and occlusion caused by the bending or twisting of the catheter. This catheter is designed to deliver CSF to the peritoneal cavity. The tubing is open-ended and contains a barium sulfate impregnated stripe to make it radiopaque throughout its entire length.

Each component of a shunting system affects the overall flow resistance of that system. To determine the flow characteristics of the enclosed shunt system, refer to the System Performance Chart.


Surgical Procedure

The introduction of a shunting system, including placement of the flushing valve, may be accomplished through a variety of surgical techniques; therefore, the surgeon is best advised to use the method which his/her own practice and training dictate to be best for the patient.

Injection

Fluids may be selectively flushed in the proximal or distal direction by injecting

Valve Performance Chart			
	Low Pressure	Medium Pressure	High Pressure
Dot Code	●	● ●	● ● ●
Closing Pressure Range	mm H ₂ O 5-50	mm H ₂ O 51-110	mm H ₂ O 111-180

Note: The above curves describe typical in-vitro valve performance characteristics for each closing pressure range. Due to characteristics of silicone materials, some variation in pressure performance may occur which historically has not compromised effective control and treatment of hydrocephalus.


Figure 4. Injection for Distal Flushing


Figure 5. Injection for Proximal Flushing

into the center of the reservoir dome and compressing the appropriate occluder. Use a 25-gauge or smaller needle for injection. **Never inject into the proximal or distal occluders or the Physiological Flow Device.**

To flush distally, introduce a 25-gauge or smaller needle into the reservoir chamber. Compressing the proximal occluder with transcutaneous finger pressure and flushing fluid into the reservoir will cause the fluid to move distally. (Figure 4.) To flush proximally, introduce a 25-gauge or smaller needle into the reservoir chamber and with transcutaneous finger pressure, compress the distal Physiological Flow Device. (Figure 5.)

A polypropylene needle-puncture shield, located in the interior of the reservoir chamber, is designed to prevent the inadvertent needle penetration of the valve base.

Patency Check

To check the patency of the shunt depress the proximal occluder with finger pressure and press down the reservoir chamber using transcutaneous finger pressure. If the reservoir chamber is not easily depressed, there may be obstruction distally. If upon releasing occluder and dome, the dome does not return easily to its original shape, there may be obstruction proximally. (Figure 6).

Pumping

For distal pumping, pump the reservoir chamber several times using transcutaneous finger pressure. (Figure 7.) For proximal pumping, compress the distal occluder (Physiological Flow Device) with finger pressure and then pump the reservoir chamber. (Figure 8.)

If the reservoir dome does not pump easily, the distal catheter may be obstructed.

If the reservoir dome does not refill promptly after pumping, the ventricular catheter or valve inlet tubing may be obstructed.

Recommendations for Use of Low-Pressure Distal Catheter. In order to allow the valve to control shunt pressure, Integra NeuroSciences recommends the use of the Integra


Figure 6. Patency Check


Figure 7. Technique for Distal Pumping


Figure 8. Technique for Proximal Pumping

NeuroSciences Peritoneal Reflux Control Catheter (Catalog NL850-1375) or Peritoneal Open-Ended Catheter (Catalog No. NL850-1376) with all pressure ranges of the **NOVUS** and **NOVUS Mini** Valves.

How Supplied

The NOVUS and NOVUS Mini Valves are supplied in a sterile, pyrogen-free, double-wrap packaging system. The double-wrap system facilitates the preferred method of sterile-product transfer from the circulating area to the sterile field. A radiopaque ventricular catheter guide is provided to help direct the catheter out of the trephine opening to minimize kinking.

Do Not Resterilize

This product is for **Single Use Only**.

Warnings

Hydrocephalic patients with cerebrospinal fluid drainage systems must be kept under close observation for signs and symptoms of increasing intracranial pressure due to shunt failure. These signs and symptoms may vary from patient to patient. Increasing intracranial pressure is characterized by headache, vomiting, irritability, listlessness, drowsiness, other signs of deterioration of consciousness, and nuchal rigidity. In the infant, increased scalp tension at the anterior fontanelle and congestion of scalp veins will be noted.

This device is made of silicone rubber, which like most rubber, may stick to itself when dry. This tendency to stick may result in some valve performance variations which may differ from label specifications. When wet, the tendency for silicone rubber to stick is reduced; therefore, the surgeon must verify that the valve components are wet and that fluid flows freely through the valve. (See **Instructions For Use**.)

This device has not been tested for drug compatibility and therefore is not intended for drug administration.

Integra NeuroSciences makes no claim for or representation as to the performance characteristics of this product if it is used in conjunction with components of other manufacturers.

Silicone tubing may be easily cut or torn when instruments are used to secure it to a connector. The use of instruments to attach silicone catheters to connectors should be avoided. When instruments are used, carefully inspect the tubing for nicks or other damage prior to closure.

Precautions

Prior to surgery, prospective patients or their representatives should be informed of the possible complications associated with the use of this product.

Fluid flow through the flushing valve should be verified immediately prior to implantation.

If a hypodermic injection into the flushing valve is required, use a 25-gauge or smaller needle and **inject through the reservoir dome only**.

The silicone proximal and distal catheter tubing should be carefully secured to the connectors with ligatures in such a manner as to avoid cutting or occluding the tubing.

Complications

Complications which may result from the use of this product include the risks associated with the medication and methods utilized in the surgical procedure, as well as the patient's response, reaction or degree of intolerance to any foreign object implanted in the body.

The principal complications associated with cerebrospinal fluid shunting into the right atrium or peritoneum are shunt obstruction, functional failure of the shunt system, infection or intracranial hypotension.

Shunt obstruction may occur in either the proximal ventricular catheter or in the distal, atrial or peritoneal catheters. Ventricular catheters may be obstructed by particulate matter such as blood clots, fibrin or brain fragments. If not properly located in the lateral ventricle, the catheter may become embedded in the ventricular wall or choroid plexus. Less commonly, the catheter may be obstructed by excessive reduction of ventricular size to slit-like proportions.

Cardiac and peritoneal catheters may

also be obstructed by particulate matter. The intra-atrial segment of the cardiac catheter may be obstructed by investment in a thrombus. Emboli from the latter may seed the pulmonary circulation sufficiently to result in pulmonary artery hypertension and cor pulmonale. Peritoneal catheters may be obstructed by the omentum or coiled loops of bowel.

Loss of valve and/or reservoir patency may result from obstruction of the fluid pathway by particulate matter such as blood clots or other biological accumulations.

Functional failure of the shunt system due to separation of its component parts can result in serious complications. Ventricular catheters may migrate into the lateral ventricles. Should the cardiac catheter become detached, it may lodge in the right atrium or ventricle or, rarely, in the pulmonary artery. Peritoneal catheters may migrate completely into the peritoneal cavity. Volvulus and perforation of intra-abdominal viscera may occur or the catheter may be extruded.

Infection is a common and serious complication of a shunting system and is most frequently caused by skin contaminants. Septicemia, which occurs most frequently in debilitated infants, can result from infections anywhere in the body and may develop with few or no symptoms. It may occur as a result of a wound infection. The presence of a foreign body (i.e. the shunting system) may trigger ventriculitis or a dormant meningitis. Intracranial infection may then be disseminated throughout the body via the distal catheter. Lesions developing from the breakdown of skin or tissue over the shunting system may also serve as foci of serious infections. In the event of an infection, removal of the shunt system is indicated in addition to the appropriate therapy.

Excessive lowering of intracranial pressure may result in complications, particularly in the infant. These include subdural hematomas, markedly sunken fontanelles, overriding of cranial bones and the conversion of a communicating to a noncommunicating hydrocephalus due to obstruction of the aqueduct of Sylvius.

Failure of the shunting system may be evidenced by any or all of the following: continuing symptoms of increased intracranial pressure, the subcutaneous exudation of CSF along the pathway of the shunt and leakage of fluid through the surgical wound. These failures require immediate replacement of the shunting system or of the affected component.

Product Information Disclosure

Integra NeuroSciences has exercised reasonable care in the choice of materials and manufacture of this product. Integra NeuroSciences excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to any implied warranties of merchantability or fitness for a particular purpose. Integra NeuroSciences shall not be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from use of this product. Integra NeuroSciences neither assumes or authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device.

Product Order Information

All products can be ordered through your Integra NeuroSciences Neuro Specialist or customer service representative or by contacting :

Integra NeuroSciences
311 Enterprise Drive
Plainsboro, NJ 08536 USA
Telephone: 1-800-654-2873
Outside the US: 1-609-275-0500
Fax: 609-275-5363
or

Integra NeuroSciences
Newbury Road, Andover
Hampshire SP10 4DR England
Tel: +44 (0) 1264-345-700
Fax: +44 (0) 1264-332-113

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. Do not use if the package has been opened or damaged.

Symbols Used On Labeling

	See instructions for use
	Expiration date
	Do not reuse after opening
	Lot number
	Sterile unless package is opened or damaged. Method of sterilization-ethylene oxide.
	Product complies with requirements of directive 93/42/EEC for medical devices.

Bibliography

Chapman PH, Cosman ER, Arnold MA. "The Relationship Between Ventricular Fluid Pressure and Body Position in Normal Subjects and Subjects with Shunts; a Telemetric Study. Neurosurgery, 26(1990), 181-189.

Foltz, Eldon L. and David B. Shurtleff. "Conversions of Communicating Hydrocephalus to Stenosis or Occlusion of the Aqueduct During Ventricular Shunt." Journal of Neurosurgery, 24 (1966), 520-529.

Hassan M, Hagashi S, Yamashita J. "Siphon-Reducing Devices in Normal-Pressure Hydrocephalus" Journal of Neurosurgery, 84 (1996), 634-641.

Ignelzi, Ronald J. and Wolff M. Kirsch. "Follow-up Analysis of Ventriculoperitoneal and Ventriculoatrial Shunts for Hydrocephalus," Journal of Neurosurgery, 42 (June, 1975), 679-682.

Illingworth, Robin D., Valentine Logue, Lindsay Symon, et al. "The Ventriculo-caval Shunt in the Treatment of Adult Hydrocephalus: Results and Complications in 101 Patients," Journal of Neurosurgery, 35 (December, 1971), 681-685.

Kuwamura, K.; Kokunal, T. "Intraventricular Hematoma Secondary to a Ventriculoperitoneal Shunt," Neurosurgery. 1982 March; 10(3); 384-386.

Little, John R., Albert L. Rhoton, Jr., and James F. Mellinger. "Comparison of Ventriculoperitoneal and Ventriculoatrial Shunts for Hydrocephalus in Children," Mayo Clin. Proc., 47 (June, 1972), 396-401.

McCullough, David C, John L. Fox, et al. "Effects of CSF Shunts on Intracranial Pressure and CSF Dynamics," Cisternography and Hydrocephalus, edited by John C. Harbert. Springfield, Illinois 1972.

Milhorat, Thomas H. Hydrocephalus and the Cerebrospinal Fluid. The Williams and Wilkins Co. Baltimore, 1972.

Portnoy, Harold D., Rudolf R. Schulte, John L. Fox, et al. "Anti-Siphon and Reversible Occlusion Valves for Shunting in Hydrocephalus and preventing Post-Shunt Subdural Hematomas," Journal of Neurosurgery, 38 (June, 1973), 729-738.

Portnoy, Harold D. and Paul D. Croissant. "Combined Drainage of Ventricular and Subdural Fluid," Surgical Neurology, 2 (January, 1974), 41-42.


Pudenz, Robert H. "The Ventriculo-Atrial Shunt," Journal of Neurosurgery, 25 (1966), 602-608.

Sugar, Oscar and Orville T. Bailey. "Subcutaneous Reaction to Silicone in Ventriculoperitoneal Shunts," Journal of Neurosurgery, 41 (September, 1974), 367-371.


Dimensional Illustrations (All dimensions are nominal)

Catalog Number	Description	Pressure
NL850-9010	Novus Valve	Low
NL850-9020	Novus Valve	Medium
NL850-9030	Novus Valve	High
NL850-9011	Novus Mini Valve	Low
NL850-9021	Novus Mini Valve	Medium
NL850-9031	Novus Mini Valve	High
NL850-9610	Novus Shunt Kit	Low
NL850-9620	Novus Shunt Kit	Medium
NL850-9630	Novus Shunt Kit	High
NL850-9311	Novus Mini Shunt Kit	Low
NL850-9321	Novus Mini Shunt Kit	Medium
NL850-9331	Novus Mini Shunt Kit	High

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© Copyright 2002 Integra NeuroSciences. All rights reserved.